



PTM450 Photocolorimetric titration LED module w/out power

Ürün Numarası:

X61T005

TRY Fiyat:
Mevcut

Bize Ulaşın

Photometric Titration Module, w/o power supply for use with OPT300 Photocolorimetric probe

Photocolorimetric titrations become fast and simple when you add a PTM450 Photocolorimetric Titration Module and OPT300 Photocolorimetric Probe to your TitrLab Titration Workstation.

Make light work of your titrations

Ideal for virtually all colorimetric titrations

Maintenance-free

Easy to handle as any standard electrode

Teknik özellikler

Ağırlık:	250 g
Boyut (Y x G x D):	45 mm x 60 mm x 109 mm
Uzunluk:	110 mm