



## C501X Cable SubD-9f/2m/SubD-9f

Ürün Numarası:

A95X501

**ÜRETİMDEN KALDIRILAN ÜRÜN**

Bu ürün artık mevcut değil.

RS232 Serial Connection Cable, SubD-9

The CCTO is a pyrex glass body cell designed for testing corrosion on large flat samples. The shape of the cell combined to the large area of the working electrode and the stainless steel counter electrode make it particularly suitable for testing low conductivity media and/or highly resistive samples. The setup includes a calomel reference electrode mounted on a vertically adjustable Luggin capillary designed to minimize ohmic drop.