



Nitrate standard, 25 mg/L NO₃⁻ (5.65 mg/L NO₃⁻-N), 500 mL

Ürün Numarası: LCW828
TRY Fiyat: Bize Ulaşın
2 hafta içinde teslimat

Standard solution for online nitrate determination

Teknik özellikler

Cihaz:	Nitratax plus sc, Nitratax clear sc, NT3100sc, NT3200sc
Hacim / Paket boyutu:	500 mL
Konsantrasyon:	25 mg/L (5,65 mg/L NO ₃ ⁻ -N)
Parametre:	Nitrate
Platform:	Çözelti

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)