



Iron Reagent Powder Pillows, 0.02-3.00mg/L Fe, FerroVer

Ürün Numarası:

85428

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Powder reagent in sealed package, for 25mL sample.

FerroVer Iron Reagent: For the determination of total soluble Iron including all ferric (Fe_3^+) and ferrous (Fe_2^+) Iron as well as complexed Iron compounds (Fe-EDTA). Perfect for many applications and a wide range of concentrations. Determines Iron in the presence of Copper.

Teknik özellikler

EPA uygunluğu:	Evet
Hacim / Paket boyutu:	1000 /pk
Metot:	FerroVer
Metot numarası:	8008
Ölçüm aralığı:	0,02 - 3,00 mg/L Fe
Platform:	Powder Pillows
Test sayısı:	1000

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)