



CuVer 2 Copper Reagent Powder Pillows

Ürün Numarası:

2188299

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For the determination of copper, 0.04-5.00mg/L Cu, in the presence of high levels of Aluminum, Calcium, and Iron. For 25 mL sample. Pack of 100.

Teknik özellikler

Hacim / Paket boyutu:	100 /pk
Metot:	Bichinchoninate
Parametre:	Copper
Platform:	Powder Pillows
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)