



Chloride Reagent Set, 0.1-25.0 mg/L Cl⁻

Ürün Numarası:

2319800

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat

Tehlikeli



Reagent set for the determination of Chloride, Mercuric Thiocyanate method.

Teknik özellikler

Cihaz:	All except DR 900
	DR/800 Series, PC II
Dipnot:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Metot:	Mercuric Thiocyanate
Metot numarası:	8113
Ölçüm aralığı:	0 - 25,0 mg/L Cl ⁻
Parametre:	Chloride
Platform:	Çözelti
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)