



Iodate-Iodide Digital Titrator Cartridge, 0.3998 N

Ürün Numarası:

1496101

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Ready-to-use titrant solution for the analysis of Sulfite, 4 - >400mg/L SO₃.

Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot:	Iodate-Iodide
Ölçüm aralığı:	0 - 800 mg/L SO ₃ ²⁻
Parametre:	Sulfite
Platform:	Dijital Titratör
Titrant:	0.3998 N KIO ₃ -KI

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)