



5500sc Reagent 3 Low Range Phosphate Package

Ürün Numarası:

6775702

TRY Fiyat:

Bize Ulaşın

Reagent 3 package for the 5500sc Low Range Phosphate analyzer. Includes 2L bottle of DI water and powder to be mixed together.

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)