



## Maintenance kit for 5500 sc Silica Analyser, 2/4 channel

Ürün Numarası:

6788302

TRY Fiyat:

Bize Ulaşın

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

### Teknik özellikler

Kutudakiler:

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Ürün Kiti:

5500 sc Silica Analyser, 2/4 Channels

### Kutudakiler

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions