



Maintenance kit for 5500sc Silica Analyser, 6 channels

Ürün Numarası:

6788303

TRY Fiyat:

Bize Ulaşın

Maintenance kit includes: air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Teknik özellikler

Kutudakiler:

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions

Ürün Kiti:

5500 sc Silica Analyser, 6 Channels

Kutudakiler

Air filter, tubing, bottle caps with colours, colorimeter cell, stir bar, fittings, pressure transducer, check valves, capillary assembly, fan filter, labels, user instructions