



CDC866T General-purpose conductivity cell, 4-pole, Temp.

Ürün Numarası: E61M015
TRY Fiyat: Bize Ulaşın
Mevcut

4-pole conductivity electrode with MAB6 plug and epoxy body. Cell constant 1cm^{-1} . Platinised, with temperature sensor. Max. temperature 80°C, min. immersion depth 35mm.

Teknik özellikler

Conductivity Cell:	CDC866T
Elektrot:	4-pole, platinised
Garanti:	12 aylar
Önemli özellik:	Conductivity electrode with MAB6 plug and epoxy body.