



Sample cells with caps for 2100 portable turbidimeters, 10 mL, 6 pcs

Ürün Numarası:

2434706

TRY Fiyat:
Mevcut

Bize Ulaşın

Replacement sample cells for use with portable turbidimeters. Pack of six round sample cells with caps.

Teknik özellikler

Ağırlık:	0,1 kg
Dipnot:	For use with Hach 2100 Series Portable Turbidimeters
Hacim / Paket boyutu:	6 /pk
İçerik:	6 sample cells with caps
Malzeme:	Glass
Numune hacmi:	10 mL