



pH test kit, 6.5 - 8.5 pH, model 17H

Ürün Numarası: 147008
TRY Fiyat: Bize Ulaşın
Mevcut

Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/Phenol red
Model:	17H
Muhafaza Türü:	B
Ölçüm aralığı:	6,6 - 8,4 pH
Parametre:	pH, mid range
Platform:	Renk Diski
Smallest increment:	0.1 pH
Test kiti modeli:	17H
Test sayısı:	200