



Test kit, phosphorus, orthophosphate (reactive) , model PO-19

Ürün Numarası:

224800

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Method / Chemistry: Color disc / Ascorbic acid. Range (mg/L): 0-1, Smallest Increment- 0.02. Range (mg/L): 0-5, Smallest Increment- 0.1. Range (mg/L): 0-50, Smallest Increment- 1. Approx. Number of Tests: 100. Case Style: 12x7.25x6 inches Blue Polypropylene Carrying Case.

Results in less than 5 minutes for clean water samples

Utilizes simple color disc comparator for accurate measurements

Model PO-19A contains filtration apparatus for turbid samples

Teknik özellikler

Metot:	Color disc/Ascorbic acid
Model:	PO-19
Muhafaza Türü:	R
Ölçüm aralığı:	Triple Range
Ölçüm aralığı (2):	0 - 4 mg/L PO ₄
Ölçüm aralığı (3):	0 - 40 mg/L PO ₄
Parameter:	Phosphorus
Parametre:	Phosphorus, Orthophosphate (Reactive) - as PO ₄ ³⁻
Platform:	Renk Diski
Smallest increment:	0.02
	0.10
	1
Standart aralık:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄
Test kiti modeli:	PO-19
Test sayısı:	100