



Test kit, phosphate, model PO-23

Ürün Numarası:

224902

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot numarası:	Color disc/Ascorbic acid
	for clear samples
Model:	PO-23
Muhafaza Türü:	F
Ölçüm aralığı:	Dual Range
Ölçüm aralığı (2):	1 - 40 mg/L PO ₄
Parameter:	Phosphate
Parametre:	Phosphate, ortho-/meta- - As PO ₄ ³⁻
Platform:	Renk Diski
Smallest increment:	0.1
	1
Standart aralık:	0.1 - 4 mg/L PO ₄ ; 1 - 40 mg/L PO ₄
Test kiti modeli:	PO-23
Test sayısı:	100