



Iron color disc test kit, model IR-18B

Ürün Numarası:

146401

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Using the simple color disc method, this kit measures Iron (medium range - 0.2-7 mg/L) and contains reagent for 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Color disc/ 1,10 Phenanthroline
Model:	IR-18B
Muhafaza Türü:	B
Ölçüm aralığı:	High Single Range
Parameter:	Iron
Parametre:	Iron, medium range - as Fe
Platform:	Renk Diski
Smallest increment:	0.25
Standart aralık:	0.25 - 7 mg/L Fe
Test kiti modeli:	IR-18B
Test sayısı:	100
Tür:	Tek Metrik