



Iron color disc test kit, model IR-18A

Ürün Numarası:

146500

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Color disc/ 1,10 Phenanthroline
Model:	IR-18A
Muhafaza Türü:	B
Ölçüm aralığı:	Low Single Range
Parameter:	Iron
Parametre:	Iron, medium range - As Fe
Platform:	Renk Diski
Smallest increment:	0.02
Standart aralık:	0.02 - 1.00 mg/L Fe
Test kiti modeli:	IR-18A
Test sayısı:	100
Tür:	Tek Metrik