



Test kit, molybdate, model MO-2

Ürün Numarası:

1419301

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Teknik özellikler

Metot:	Color disc/ Mercaptoacetic acid
Model:	MO-2
Muhafaza Türü:	F
Ölçüm aralığı:	Dual Range
Ölçüm aralığı (2):	0 - 40 mg/L Mo
Parametre:	Molybdate - As Molybdenum
Platform:	Renk Diski
Smallest increment:	0.2
	1
Test kiti modeli:	MO-2
Test sayısı:	100