



Test kit, manganese color disc, model MN-PAN

Ürün Numarası:

2350800

Tehlikeli

TRY Fiyat:

Bize Ulaşın

2 hafta içinde teslimat



Manganese causes a bitter taste in water, and at concentrations above 0.1 mg/L, it causes objectionable stains on laundry and plumbing fixtures. Manganese control also is important in the beverage, paper, and textile industries, and in dye production and food processing. Utilizing the simple color disc method, this kit measures manganese (Range: 0.05-0.7 mg/L) and contains 50 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Metot:	Color disc/PAN
Model:	MN-PAN
Muhafaza Türü:	D
Ölçüm aralığı:	0 - 0,70 mg/L Mn
Parameter:	Manganese
Parametre:	Manganese - As Mn
Platform:	Renk Diski
Smallest increment:	0.05
Standart aralık:	0 - 0.70 mg/L Mn
Test kiti modeli:	MN-PAN
Test sayısı:	50