



## Test kit, glutaraldehyde, model GT-1

Ürün Numarası:

2587200

Tehlikeli

TRY Fiyat:  
Mevcut

Bize Ulaşın



From oil fields to cooling towers, hospitals to food plants, poultry farms to dental offices, the Glutaraldehyde Test Kit can help monitor disinfection processes. The test kit measures form 0.5 to 4000 ppm glutaraldehyde.

**Helps monitor disinfection processes**

**Simple, low-cost color disc method is the most accurate visual test kit method**

**Achieve fast results with an easy-to-carry test kit**

---

### Teknik özellikler

Metot:	Color disc
Model:	GT-1
Muhafaza Türü:	Q
Ölçüm aralığı:	1 - 2800 mg/L multiple ranges
Parameter:	Glutaraldehyde
Parametre:	Glutaraldehyde
Platform:	Renk Diski
Smallest increment:	0.5
Standart aralık:	1 - 2800 mg/L multiple ranges
Test kiti modeli:	GT-1
Test sayısı:	100