



Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Ürün Numarası:

143601

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO₂, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Teknik özellikler

Dipnot:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Metot:	Drop count titration/ Sodium hydroxide Phenolphthalein
Muhafaza Türü:	D
Ölçüm aralığı:	2 - 40 mg/L
Ölçüm aralığı (2):	1.25 - 25 mg/L
Ölçüm aralığı (3):	5 - 100 mg/L CO ₂
Parametre:	Carbon Dioxide - As CO ₂
Platform:	Drop Count Titration
Smallest increment:	1.25 2 5
Test kiti modeli:	CA-23
Test sayısı:	200