



Test kit cyanide, model CYN-3

Ürün Numarası:

201002

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat

Tehlikeli



The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures.

Simple, low-cost method using a stepped color comparator

Achieve fast results with an easy-to-carry test kit

Teknik özellikler

Metot:	Color disc/ Pyridine-Pyrazolone
Model:	CYN-3
Muhafaza Türü:	R
Ölçüm aralığı:	0 - 0,3 mg/L CN ⁻
Parameter:	Cyanide
Parametre:	Cyanide - As free CN ⁻
Platform:	Renk Diski
Smallest increment:	0.02
Standart aralık:	0 - 0.3 mg/L CN ⁻
Test kiti modeli:	CYN-3
Test sayısı:	100