



## Hardness (Total & Calcium) Test Kit, Model HAC-DT

Ürün Numarası:

2063900

TRY Fiyat:  
Mevcut

Bize Ulaşın

Tehlikeli



Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

### Teknik özellikler

Kutudakiler:	Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.
Metot:	Digital Titrator/EDTA
Muhafaza Türü:	F
Ölçüm aralığı:	10 - 4,000 mg/L $\text{CaCO}_3$
Parametre:	Hardness, total & Calcium - As $\text{CaCO}_3$
Platform:	Dijital Titratör
Smallest increment:	0.1-10
Test kiti modeli:	HAC-DT
Test sayısı:	100 Ca +100 T

### Kutudakiler

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.