



Test kit, alkalinity, model AL-TA

Ürün Numarası:

2314500

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Dipnot:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Metot:	Drop count titration/Sulfuric acid
Muhafaza Türü:	B
Ölçüm aralığı:	23 - 495 gpg & 385 - 8500 mg/L
Parametre:	Alkalinity - (P) & Total (MO) As CaCO ₃
Platform:	Drop Count Titration
Smallest increment:	23 gpg
	385
Test kiti modeli:	AL-TA
Test sayısı:	100