



Test kit, alkalinity, model AL-AP

Ürün Numarası:

2444300

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Teknik özellikler

Dipnot:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. **1 gpg = 17.1 mg/L (MO) = Methyl Orange
Metot:	Drop count titration/Sulfuric acid
Muhafaza Türü:	B
Ölçüm aralığı:	0 - 8 gpg & 1 - 20 gpg
Parametre:	Alkalinity - (P) & Total (MO) As CaCO ₃
Platform:	Drop Count Titration
Smallest increment:	0.4 gpg
	1 gpg
Test kiti modeli:	AL-AP
Test sayısı:	100