



## Test kit, triazole, model TZ-1, w/ 230 vac uv lamp & power supply

Ürün Numarası:

2167502

TRY Fiyat:  
2 hafta içinde teslimat

Bize Ulaşın

Tehlikeli



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

**Measures common corrosion and biological growth inhibitor using simple color disc method**

**UV lamp digestion**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

### Teknik özellikler

Güç kaynağı:	230 Vac UV Lamp and Power Supply
Metot:	Color disc/ UV photolysis (230 Vac, 50 Hz)
Model:	TZ-1
Muhafaza Türü:	F
Ölçüm aralığı:	0 - 15,0 mg/L
Parametre:	Triazole -As benzotriazole
Platform:	Renk Diski
Smallest increment:	0.5
Test kiti modeli:	TZ-1
Test sayısı:	50