



Bottle, BOD, 300 mL, 24/pk

Ürün Numarası: 62124
TRY Fiyat: Bize Ulaşın
Mevcut

The standard for BOD analysis. Borosilicate clear glass bottles with flared lip and ground-glass stoppers to prevent air entrapment and for water seal. Frosted marking area. Stoppers included. Pack of 24.

Teknik özellikler

Body Material:	Glass
Hacim / Paket boyutu:	24 /pk
Kapasite:	300 mL
Malzeme:	Borosilicate Glass
Package size UOM:	24pk