



Beaker, polypropylene, low form, 400 mL, 6/pk

Ürün Numarası: 108078
TRY Fiyat: Bize Ulaşın
Mevcut

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Teknik özellikler

Birim:	6 /pk
Çap:	99.1 mm with spout
Gövde malzemesi:	Polypropylene
Kapasite:	400 mL
Malzeme:	Polypropylene, Low Form