



CEL/890 Advanced portable laboratory

Ürün Numarası:

2994000

Tehlikeli

ÜRETİMDEN KALDIRILAN ÜRÜN
Bu ürün artık mevcut değil.



Offering on-site results in an easy-to-use, rugged, portable lab, Hach's Oil & Gas Water Quality Laboratory covers parameters critical to effective source water and base water analysis. Depending on water composition, this kit may also be used to support analysis of produced water, flowback water, water treatment, and enhanced oil recovery applications.

Teknik özellikler

Kutudakiler:

Kit contains a DR900 Colorimeter, HQ40D Multi-Parameter Digital Meter, Digital Titrator, all required reagents for 12 parameters critical to oil and gas measurement, apparatus and accessories required to perform tests, and two convenient, rugged carrying cases.

Uygulama:

Fracing

Kutudakiler

Kit contains a DR900 Colorimeter, HQ40D Multi-Parameter Digital Meter, Digital Titrator, all required reagents for 12 parameters critical to oil and gas measurement, apparatus and accessories required to perform tests, and two convenient, rugged carrying cases.