



Be Right™



Rectangular sample cell, 1"x1", 25 mL, glass (2 pieces)

Ürün Numarası:

2095000

TRY Fiyat:

Bize Ulaşın

Mevcut

Accessories for spectrophotometers/spectral measurements of liquids.

For use with the DR2010 and DR2000 spectrophotometers. Comes as a set of matched pair cells.

Teknik özellikler

Hacim / Paket boyutu:

2

Malzeme:

Glass

Optik yol uzunluğu:

1 inch rectangular, 25 mL