



Nitrate standard, 100 mg/L NO₃ (22.6 mg/L NO₃-N), 500 mL

Ürün Numarası:

LCW826

TRY Fiyat:
Mevcut

Bize Ulaşın

Standard solution for online nitrate determination

Teknik özellikler

Cihaz:	Nitratax plus sc, Nitratax clear sc, NT3100sc, NT3200sc
Hacim / Paket boyutu:	500 mL
Konsantrasyon:	100 mg/L NO ₃ (22.6 mg/L NO ₃ -N)
Parametre:	Nitrate
Platform:	Çözelti

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)