



## MolyVer 3 Molybdenum Powder Pillows (3), pk/100

Ürün Numarası:

1417869

Tehlikeli

TRY Fiyat:  
Mevcut

Bize Ulaşın



Powder reagent in sealed package for the determination of high range molybdenum (0-40.0mg/L Mo), Mercaptoacetic Acid Method. For 25mL sample.

### Teknik özellikler

|                       |                     |
|-----------------------|---------------------|
| Hacim / Paket boyutu: | 100 /pk             |
| Metot:                | Mercaptoacetic Acid |
| Metot numarası:       | 8036                |
| Parametre:            | Parameter           |
| Platform:             | Powder Pillows      |
| Test sayısı:          | 100                 |

### Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)