



Silver nitrate digital titrator cartridge, 0.2256 N, 1 piece

Ürün Numarası:

1439601

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Ready-to-use titrant solution for the analysis of Chloride (Silver Nitrate Method), 10 - 100 mg/L Cl⁻.

Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot:	Silver Nitrate
Metot numarası:	8207
Ölçüm aralığı:	10 - 100 mg/L Cl ⁻
Parametre:	Chloride
Platform:	Dijital Titratör
Test sayısı:	40 - 100
Titrant:	0.2256 N AgNO ₃

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)