



Sodium Thiosulfate Digital Titrator Cartridge, 0.113 N

Ürün Numarası:	2267301
TRY Fiyat:	Bize Ulaşın
Mevcut	

Ready-to-use titrant solution for the analysis of Total Chlorine, 20-2000 mg/l Cl₂.

For chlorine determination by iodometric Digital Titrator titration. 0.113 N, Digits x 0.1 = mg/L Cl₂ in a 5 mL sample, Digits x 0.000141 = meq. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot:	Iodometric
Metot numarası:	8209
Ölçüm aralığı:	20 - 2000 mg/L Cl ₂
Parametre:	Chlorine, Total
Platform:	Dijital Titratör
Test sayısı:	40 - 100
Titrant:	0.113 N Na ₂ S ₂ O ₃

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)