



Reagent powder pillows, DPD total chlorine, 5 mL, pk/1000

Ürün Numarası:

1407628

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For determination of total chlorine by the DPD method.

For use with Hach's Chlorine Visual Test Kits.

Sample Size: 5 mL. Range: up to 3.5 mg/L Cl_2 .

Pack of 1000 powder pillows.

Boric Acid free!

Teknik özellikler

Metot:	DPD
Numune boyutu:	5 mL
Platform:	Powder Pillows
Saklama koşulları:	10 - 25 °C
Test sayısı:	1000

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)