



Free Chlorine Reagent Powder Pillows, DPD (pk/100)

Ürün Numarası:

1407799

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For determination of free chlorine by the DPD method. For use with Hach's Chlorine Visual Test Kits.

Sample Size: 5 mL. Range: up to 3.5 mg/L Cl₂.

Pack of 100 powder pillows.

Boric Acid free!

Teknik özellikler

Metot:	DPD
Numune boyutu:	5 mL
Parametre:	Chlorine, Free
Platform:	Powder Pillows
Saklama koşulları:	10 - 25 °C
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)