



DPD Free Chlorine Reagent Powder Pillows, 5 mL, pk/50

Ürün Numarası:

2197846

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For determination of free chlorine with Hach's Chlorine Visual Test Kits- 2060300. Sample Size: 5 mL. Range: up to 2.5 mg/L Cl_2 . Pack of 50 Powder Pillows.

Boric Acid free!

Teknik özellikler

Metot:	DPD
Numune boyutu:	5 mL
Parametre:	Chlorine, Free
Platform:	Powder Pillows
Saklama koşulları:	10 - 25 °C
Test sayısı:	50

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)