



Reagent set, sulfite, digital titrator

Ürün Numarası:

2272300

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



For determination of sulfite by iodate-iodide Digital Titrator titration, Hach Method 8216. Range: 4 - 400 mg/L SO₃. Reagent Set includes: Dissolved Oxygen 3 Reagent Powder Pillows, Iodate-Iodide Cartridge (0.3998 N), and Starch Indicator. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot:	Iodate-Iodide
Metot numarası:	8216
Önemli özellik:	Sulfite Reagent set, Digital Titrator
Platform:	Dijital Titrator
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)