



Calibration set, Stablcal turbidity standards, 2100P portable turbidimeter, 100 mL bottles

Ürün Numarası:

2659410

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Four ready to use, stabilised Formazine Primary standard solutions in bottles. Kit includes 100 mL of each: <0.1 NTU, 20 NTU, 100 NTU, 800 NTU .

Teknik özellikler

Hacim / Paket boyutu:

Set of 4

Konsantrasyon:

< 0,1, 20, 100, 800 NTU

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)