



NitraVer 5 Nitrate Reagent PP, 5 mL pk/100

Ürün Numarası:

1403599

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Powder reagent in sealed packages for the determination of Nitrogen-Nitrate up to 50 mg/l NO₃-N.

For high range nitrate determination by the NitraVer 5 cadmium reduction method. Range: up to 50 mg/L as nitrate-N. Sample Size: 5 mL. Pack of 100 powder pillows.

Teknik özellikler

Hacim / Paket boyutu:	100 /pk
Metot:	Cadmium Reduction
Platform:	Powder Pillows

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)