



SPADNS 2 Fluoride Reagent, arsenic-free, ACCUVAC, 0.02-2.00 mg/L F

Ürün Numarası:

2527025

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



ACCUVAC ampules contain the precise amount of reagent for a single test and can be used as a measurement cuvette. For the determination of Fluoride by the SPADNS 2 method. pk/25

Teknik özellikler

EPA uygunluğu:	Evet
Hacim / Paket boyutu:	25 /pk
Method Number 1 Chemistry:	10225
Metot:	SPADNS 2
Ölçüm aralığı:	0,02 - 2,00 mg/L F ⁻
Parametre:	Fluoride
Platform:	Accuvac
Test sayısı:	25

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)