



Citric Acid / Surfactant Reagent for Silica analysis, 2.9 L

Ürün Numarası:

2347003

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Citric Acid / Surfactant Reagent (2.9 L) used in the Series 5000 Silica Analyser.

Teknik özellikler

Cihaz:	5500sc Silica
Hacim / Paket boyutu:	2.9 L
Ölçüm aralığı:	N/A
Parametre:	Silica
Saklama koşulları:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)