



Amino Acid F Reagent for Silica analysis, 2.9 L

Ürün Numarası:

2353103

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Amino Acid F Reagent (2.7 L) used in the Series 5000 Silica Analyser.

Amino Acid F reagent is supplied in two parts to ensure a long shelf life: a 3-liter bottle of Amino Acid F Solution and a bottle of Amino Acid F Powder. Powder must be added to the solution immediately before installation in the analyser.

Teknik özellikler

Cihaz:	5500sc Silica
Hacim / Paket boyutu:	1 x 2.7 L & 1 x 308 g
Ölçüm aralığı:	N/A
Parametre:	Silica
Saklama koşulları:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)