




Stablcal Turbidity Standard Bulk Solution Calibration Kit, 20 NTU and 0.1 NTU

Ürün Numarası:

2659600

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Stablcal Turbidity Standard Bulk Solution Calibration Kits for Hach 1720 Series Turbidimeters. USEPA accepted for calibration purposes. Stablcal Stabilised Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. Stablcal Kit includes 4 each, 1 liter bottles of: 0.1 NTU, 20 NTU.

Precise, Pre-mixed Formazin Solutions

Stablcal Stabilised Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations as low as 0.060 NTU for conventional turbidimeters. Shelf life is guaranteed for two years from the date of manufacture. Whether you are monitoring wastewater effluent, high-quality treated drinking water or an industrial process stream, you can be completely confident in the results you obtain with Hach Stablcal Stabilised Formazin Standards.

No special glassware or difficult lab techniques

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Since then, Formazin has been the only primary standard that is accepted virtually throughout the world for regulatory reporting. All other turbidity standards, both alternate and secondary, must be traced to it.

Teknik özellikler

Hacim / Paket boyutu:	4x 1 L each
Konsantrasyon:	20 NTU
	0.1 NTU
Ölçüm aralığı:	N/A
Parametre:	Turbidity
Saklama koşulları:	Store between 5 - 25 °C, away from direct sunlight

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)