



Hardness (Calcium and Total) Reagent Set, Digital Titrator

Ürün Numarası:

2272100

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



For Calcium and Total Hardness determination by EDTA Digital Titrator titration. Ranges: 10-160 and 100-4000 mg/L CaCO₃. Reagent Set contains: Buffer Solution (Hardness 1), CalVer 2 Calcium Indicator Powder Pillows, ManVer 2 Hardness Indicator Powder Pillows, Potassium Hydroxide Solution (8 N), EDTA Titration Cartridges (1 each) (0.800 and 0.0800 M), and 5.25N sulfuric acid.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot:	EDTA Titration
Ölçüm aralığı:	10 - 160 mg/L & 100-4000 mg/L CaCO ₃
Önemli özellik:	Reagent set, Digital Titrator
Parametre:	Hardness (Calcium and Total)

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)