



Filter paper, glass fiber, extra thick, pore size 1 µm, diameter 47 mm, 100/pk

Ürün Numarası:

2842100

TRY Fiyat:

Bize Ulaşın

Mevcut

For filtration of heavily contaminated samples. High particulate-holding capacity makes discs efficient depth filters and allows for filtration of large volumes of solutions. Pall Life Sciences.

Teknik özellikler

Çap:	47 mm mm
Gözenek boyutu:	1 µm
Hacim / Paket boyutu:	100
Malzeme:	Glass Microfiber