



3M KCl electrolyte for reference electrode, 500mL

Ürün Numarası: Z363140,00500
TRY Fiyat: Bize Ulaşın
Teslimat tarihi için arayın

Consumables for a sodium analyser in its spare kit.

Teknik özellikler

Hacim / Paket boyutu: 500 mL
Konsantrasyon: 3.0 mol/L
Önemli notlar: Reference electrolyte for refillable reference electrodes

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)