



Bottle BOD, 300 mL, disposable

Ürün Numarası:

2943101

TRY Fiyat:
Mevcut

Bize Ulaşın

Single 300 mL disposable BOD bottle without stopper. Plastic bottle is EPA approved for BOD measurement.

Box of 100 bottles available under article number 2943100. Plastic stoppers (25 pack) available under article number 2943900. Glass stoppers (quantity of 1) available under article number 190901.

Teknik özellikler

Body Material:	Plastic
Hacim / Paket boyutu:	1 each
Kapasite:	300 mL
Malzeme:	PET