



Reagent set, hardness (total), HA-71A test kit

Ürün Numarası:

2437800

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



For determination of Total Hardness by the Drop Count EDTA Titration Method. Replacement reagent set for HA-71A Test Kit. Set contains Buffer Solution (Hardness 1), Hardness 2 and Hardness 3 Test Solutions for approximately 100 tests.

Teknik özellikler

Metot:	EDTA Titration
Ölçüm aralığı:	1 - 20 mg/L CaCO ₃
Önemli özellik:	Reagent set, HA-71A test kit
Parametre:	Hardness (Total)

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)