



Hardness (Calcium & Total Sequential) LR Reagent Set, Digital Titrator

Ürün Numarası:

2448600

Tehlikeli

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat



For sequential determination of Calcium and Total Hardness by EDTA Digital Titrator titration. Hach method 8329. Range: 10-160 mg/L CaCO₃. Set contains Buffer Solution (Hardness 1), CalVer 2 Calcium Indicator Powder Pillows, ManVer 2 Hardness Indicator Powder Pillows, Potassium Hydroxide Solution (8 N), Sulfuric Acid (5.25 N) and 0.0800 M EDTA Digital Titrator Cartridge. Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot:	EDTA Titration
Metot numarası:	8329
Ölçüm aralığı:	10 - 160 mg/L CaCO ₃
Önemli özellik:	Reagent set, LR, Digital Titrator
Parametre:	Hardness (Calcium and Total Sequential)
Platform:	Dijital Titratör
Test sayısı:	40 - 100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)